

JUL 22, 2026

FROM BENCHMARKING TO AGENTIC WORKFLOWS

Foundation Models For X-ray Data Analysis

XIANGYU YIN

Assistant Computational scientist
Computational Science and AI Group
Advanced Photon Source
Argonne National Laboratory




Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.




**Advanced
Photon Source**

Ptychography



Retrieve object $O(\vec{r})$ (and phase $P(\vec{r})$) from a series of scanning diffraction patterns $(|\tilde{\psi}(\vec{r}_p, \vec{k})|^2)$

Sensitive to reconstruction parameters and prone to artifacts



Sample: integrated circuit


Instrument: 2x_fm

Diagnosis:

- Local scan drift

Actions:

- Turn on position correction
- Crop data to 64x64
- Increase the number of iterations



Sample: Ecoli


Instrument: Bionanoprobe

Diagnosis:

- Slowly-varying background
- Noisy artifacts

Actions:

- Use variable probe correction
- Use Poisson noise model



Sample: super crystal

Instrument: Ptycho-SAXS

Diagnosis:

- Non-uniform background
- Singularity artifacts in the object

Actions:

- Use compact batch scheme with momentum acceleration
- Multislice ptychography

Benchmark LLMs/VLMs on Ptychography Tasks

Task 1: Visual artifact detection

Find artifacts (e.g., grid/aliasing/ring) in reconstructed ptychography images.

Input: reconstructed image

Output: diagnosis

Task 2: Parameter recommendation

Given a dataset's metadata, recommend new/better reconstruction parameters.

Input: sample metadata


Output: actions

Strategies compared

In-context learning (ICL): prompting w. examples

Sample-specific few-shot (SSFS): fetch relevant examples

Supervised fine-tuning (SFT): train adaptors w. few examples



Sample: integrated circuit


Instrument: 2x_fm

Diagnosis:

- Local scan drift

Actions:

- Turn on position correction
- Crop data to 64x64
- Increase the number of iterations



Sample: Ecoli


Instrument: Bionanoprobe

Diagnosis :

- Slowly-varying background
- Noisy artifacts

Actions:

- Use variable probe correction
- Use Poisson noise model



Sample: super crystal

Instrument: Ptycho-SAXS

Diagnosis :

- Non-uniform background
- Singularity artifacts in the object

Actions:


- Use compact batch scheme with momentum acceleration
- Multislice ptychography



U.S. DEPARTMENT
of ENERGY

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.


Benchmark LLMs/VLMs on Ptychography Tasks



Task 1: Visual artifact detection

Best recipe:
SFT + sample-specific few-shot examples

Example relevancy matters most,
irrelevant context can hurt performance




Task 2: Parameter recommendation

Best recipe:
large base + sample-specific few-shots

SFT can produce “rigid expert”: great 0-shot but few shots degrades performance

Build Workflows around Tested Models and Recipes




Case Study: AI-Assisted 2D Materials Recon.

arXiv:2410.09034

AI workflow can improve
reconstruction quality in
a few iterations

Bayesian optimization
can take tens/hundreds
of iterations



Human-in-the-loop Result Feedback


PEAR: Your command execution is successful :) Please check your reconstructions and answer some questions:
PEAR: Did the reconstruction converge?
User: yes
PEAR: Do you see any grid artifacts in the object?
User: no
PEAR: Is the initial probe accurate enough?
User: yes
PEAR: Do you see any structures in the last probe mode?
User: no
PEAR: Is there anything else you want to tell me?
User: each layer seems to have random features
PEAR: Thank you for your feedback!

Parameter Recommendation

PEAR: Summarizing previous conversation...
The reconstruction converged successfully.

- There are no grid artifacts present in the object.
- The initial probe is accurate enough.
- There are no visible structures in the last probe mode.
- Each reconstructed layer has random features

Recommendation:
* Increase the layer regularization by 0.3.



Human-in-the-loop Result Feedback


PEAR: Your command execution is successful :) Please check your reconstructions and answer some questions:
PEAR: Did the reconstruction converge?
User: yes
PEAR: Do you see any grid artifacts in the object?
User: no
PEAR: Is the initial probe accurate enough?
User: yes
PEAR: Do you see any structures in the last probe mode?
User: no
PEAR: Is there anything else you want to tell me?
User: the atoms still look a bit blurred
PEAR: Thank you for your feedback!

Parameter Recommendation

PEAR: Summarizing previous conversation...
The reconstruction converged successfully.

- There are no grid artifacts present in the object.
- The initial probe is accurate enough.
- There are no visible structures in the last probe mode.
- The atoms still appear a bit blurred.

Recommendation:
** Increase the standard deviation of the gaussian kernel for detector blurring by 0.5.



Human-in-the-loop Result Feedback

PEAR: Your command execution is successful :) Please check your reconstructions and answer some questions:
PEAR: Did the reconstruction converge?
User: yes
PEAR: Do you see any grid artifacts in the object?
User: no
PEAR: Is the initial probe accurate enough?
User: yes
PEAR: Do you see any structures in the last probe mode?
User: no
PEAR: Is there anything else you want to tell me?
User: no
PEAR: Thank you for your feedback!


Parameter Recommendation

PEAR: Summarizing previous conversation...
The reconstruction converged successfully.

- There are no grid artifacts present in the object.
- The initial probe is accurate enough.
- There are no visible structures in the last probe mode.
- The user has no additional information

Recommendation:
** Increase the standard deviation of the gaussian kernel for detector blurring by 0.5.

From Parameter Tuning to Algorithm Discovery



Can you help me find a regularizer to suppress grid artifacts in IC dataset





Action: GENERATE

```
def regularize_llm(self):  
    """Physics-aware  
    regulariser for hard-X-ray  
    IC ptycho."""  
    <function content>  
    self.set_data(updated_obj)
```

Here are the reconstructed object and probe:

Code
Execution



Research
Agent

Coding
Agent

Evaluation
Agent

Analysis
Agent

Below is a compact "look-up sheet" you can keep next to your reconstruction script. It lists the regularizers, constraints and mitigation steps that have proved useful for hard-X-ray ptychography of planar integrated-circuit (IC) samples acquired in step-scan mode with (i) an inaccurate initial probe, (ii) very sparse/large scan steps and (iii) visible raster-grid artifacts.

Regularization techniques that work for this experiment

A. Spatial-domain sparsity & piece-wise constancy

- Isotropic Total-Variation (TV) or Total-Generalised-Variation (TGV) penalties suppress noise and fill in missing pixels when the overlap is small; TV was shown to recover 100 nm IC lines at only 15 % overlap (arxiv.org).
- Structure-Tensor or directional TV priors preserve the long, straight metal lines typical of IC routings (arxiv.org).

B. Frequency-domain / smoothness penalties

- L2/Tikhonov or gradient-penalising terms on the object and/or probe, implemented in modern autodiff engines (e.g. phaser/ptyRAD), keep high-frequency probe noise from leaking into the object (arxiv.org, ar5iv.labs.arxiv.org).

Action: EVOLVE




```
def regularize_llm(self):  
    """Multi-scale anisotropic  
    DTGV-Huber regulariser with  
    adaptive notch filtering."""  
    <function content>  
    self.set_data(updated_obj)
```

- **(Ground truth) metrics mode**
- **Human-in-the-loop mode**
- **Auto mode (tools & VLMs)**


- Expose the main hyper-parameters (gamma, lam_init/min/max, huber_eps_factor, phase_tv_strength, softplus_beta_scale) through the external API so they can be tuned per dataset; grid/Bayesian search typically yields ≥ 2 dB PSNR
- Keep the Barzilai-Borwein step length as a tensor (do not call .item()) and compute it in full complex arithmetic; this preserves autograd graphs for unrolled or meta-learning pipelines.

...


Discovered Regularizer v.s. the Hand-Crafted Baseline



Discovered




Baseline




Algo. Evolution History Recorded for Insights

Algorithm Discovery Evolution



Technique Timeline: When Techniques Appear Across Generations



Key Technical Insights:

Adaptive notch filtering:

FFT-based grid artifact suppression

Structure-tensor TV:

Edge-aware anisotropic smoothing

Barzilai-Borwein:


Adaptive step size acceleration

Most Data Analysis Problems are Unstructured




Ming Du

Structured problem




Challenges for unstructured analysis



Development environment setup

Evaluation script

Unstructured problem




Troubleshooting

Exploration-exploitation tradeoff


Automated Agentic System for Unstructured Analysis



Ming Du



Example: finding a translational registration algorithm for small, noisy fluorescence images



Algorithm	Development set	Holdout set
Brute-force error minimization	1.25	0.98
Phase correlation	5.28	5.59
OpenEvolve	0.20	0.23
ERA	0.46	0.50
CVEvolve	0.11	0.12



Du, M., Yin, X., Luo, Y., Beniwal, D., Tang, S., Sharma, H., & Cherukara, M. J. (2026). CVEvolve: Autonomous Algorithm Discovery for Unstructured Scientific Data Processing. *arXiv preprint arXiv:2605.11359*.



<https://github.com/AdvancedPhotonSource/CVEvolve>

Skills + General Agent for Interactive Flexible Analysis



Yi Jiang

Ptycho Skills

If you need...	See...
Data formats & required HDF5 keys	references/data_formats.md
Parameter reference & suggestions (pear wrapper params)	references/ptychox_parameters.md
Troubleshooting artifacts & failures	references/troubleshooting.md
Quality assessment & "what to change next"	references/quality_assessment.md
Physics background & key equations	references/theoretical_foundations.md
Deep API customization (Options, Schedule, constraints, losses)	references/ptychox_api_guide.md
Warm-start optimization (multi-stage refinement)	references/warm_start_guide.md
Batch reconstruction (many scans)	references/batch_reconstruction.md
X-ray template (step-scan)	templates/ptychox_xray_standard.py
X-ray template (fly-scan)	templates/ptychox_flyscan.py
Near-field template	templates/ptychox_nearfield.py
Multislice template	templates/ptychox_multislice.py
Electron (4D-STEM) template	templates/ptychox_electron_stem.py
Custom regularizer template	templates/ptychox_custom_regularizer.py
JSON schema & sample configs	examples/
Running on ALCF (Polaris/Aurora)	See "Running on ALCF" section below
What makes PtychoX different (vs pty-chi/ptyrad/phaser/foldslice)	docs/unique_features.md
PtychoX vs pty-chi (algorithm-level differences, ported bug fixes)	docs/vs_ptychi.md

Self-updating of Ptycho Skills

"Can you process a new data?"



"Explore more parameters/strategies to improve image quality"



"Summarize the lessons we learned and update the Ptycho Skill"



U.S. DEPARTMENT
of ENERGY

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.

Argonne
NATIONAL LABORATORY



Advanced
Photon Source

The AI era in X-ray data science is not about replacing physics and scientists with black boxes.

It is about rebuilding data analysis systems so physics, AI, and human expertise work together.

Thanks! I look forward to the discussion.

Xiangyu Yin

*Assistant Computational Scientist
Computational Science and AI Group
Advanced Photon Source
Argonne National Laboratory*

Acknowledgement: Arthur Glowacki, Yi Jiang, Wendy Di, Grace Luo, Olga Antipova, Luxi Li, Si Chen, Junjing Deng, Barry Lai, Juanjuan Huang, Shelly Kelly, Peco Myint, Michael Prince, Sinisa Vesali, Nicholas Schwartz, Ming Du, Aileen Luo, Nina Andrejevic, Mathew Cherukara, Stefan Vogt, Alec Sandy, Jonathan Lang, Laurent Chapon

This research used resources of the Advanced Photon Source, a U.S. Department of Energy (DOE) Office of Science User Facility operated for the DOE Office of Science by Argonne National Laboratory under Contract No. DE-AC02-06CH11357. I acknowledge funding support from Argonne **LDRD 2023-0049**